

Cybertech Nitrile Examination Gloves

Enhanced Grip and Protection.



TO ORDER VISIT:
KENTEXPRESS.CO.UK




CYBERTECH™
Core Products. Exceptional Value.

Cybertech Nitrile Examination Gloves

- Microtextured fingertips ensure a secure grip on instruments and materials
- Beaded cuff prevents rolling back, 0.07mm thickness at fingertips/0.05mm thickness at palm
- Powder free to reduce the risk of skin and respiratory sensitivities
- Soft, comfortable ambidextrous design
- Single-use
- Non-sterile
- Not made with natural rubber latex



CLASS 1 MEDICAL DEVICE CONFORMING TO REGULATION
EU MDR 2017/745 AND COMPLY TO EN 455-1/2/3/4
AND ASTM D6319 • AQL 1,5



ITEM DESCRIPTION

QUANTITY

ITEM CODE

Cybertech Nitrile Examination Gloves, XS	100/Box	(5724323)
Cybertech Nitrile Examination Gloves, S	100/Box	(5724324)
Cybertech Nitrile Examination Gloves, M	100/Box	(5724325)
Cybertech Nitrile Examination Gloves, L	100/Box	(5724326)
Cybertech Nitrile Examination Gloves, XL	100/Box	(5724327)

TO ORDER VISIT: [KENTEXPRESS.CO.UK](https://www.kentexpress.co.uk)

